

Non WHO Reference Material NIBSC Anti-CMV IgM Serology Reagent NIBSC code: 25/B919 Instructions for use (Version 2.0, Dated 24/01/2025)

This material is not for in vitro diagnostic use

1. INTENDED USE

This reagent is intended for performance monitoring of immunoassays detecting antibodies to CMV IgM.

The reagent is expected to be detected by assays with high CMV IgM assays, different batches of the same assay kit, and/or different assay conditions, may give different values with this reagent.

The reagent is not intended to be used for assay calibration purposes, or for in vitro diagnostic use.t

2. CAUTION

This preparation is not for administration to humans or animals in the human food chain.

25/B919 has been prepared from a pool of anti – CMV IgM reactive defibrinated plasma donations, repeatedly reactive in commercial EIA kits . The reactive donations used to prepare 25/B957 were non-reactive for HBsAg, anti-HIV1/2 and anti-HCV using commercial EIA kits. The reactive donations were pooled and then diluted in a pool of defibrinated human plasma donations. These samples were non-reactive for HBsAg, anti-HCV and anti-HIV1/2, Anti-HTLV and Anti-Syphilis using commercial EIA kits. Bronidox[®] was added toa concentration of 0.05%(w/v) as a preservative. As with all materials of biological origin, this preparation should be regarded as potentially hazardous to health. It should be used and discarded according to your own laboratory's safety procedures. Such safety procedures should include the wearing of protective gloves and avoiding the generation of aerosols. Care should be exercised in opening ampoules or vials, to avoid cuts.

3. UNITAGE

There is no unitage assigned to this product.

4. CONTENTS

Country of origin of biological material: United Kingdom. Ready-to-use reagent REF 25/B919 1x7mL Blood Tubes Defibrinated Plasma 1mL Bronidox[®] (Sigma-Aldrich) 0.05% (w/v)

5. STORAGE

Reagents are to be kept at 2-8°C upon receipt

o Reagents may be stored at 2-8°C until use by date

o For single use only reagents should be divided into measured aliquots of one use and stored below -20°C to avoid freeze/thaw cycles. Once thawed use immediately.Do not refreeze

o Ensure all containers are properly sealed to avoid drying out of the reagent o Avoid microbial contamination of this product as this may alter product performance

o Avoid excessively high temperatures or humidity

Material type: Liquid – will be shipped according to the storage and shipping conditions of the product



6. DIRECTIONS FOR OPENING

Vials have a screw cap; an internal stopper may also be present. The cap should be removed by turning anti-clockwise. Care should be taken to prevent loss of the contents. Please note: If a stopper is present on removal of the cap, the stopper should remain in the vial or be removed with the cap.

7. USE OF MATERIAL

. Use of this reagent is to be restricted to trained laboratory staff only

- 2. Use suitable (latex/nitrile) gloves and eye/skin protection
- 3. Include reagent as a normal sample in routine work list
- 4. Allow reagent to reach room temperature before use

5. Plot reagent result on a QC chart to monitor performance or using NIBSC Results reporting Software.

8. STABILITY

Reference materials are held at NIBSC within assured, temperature-controlled storage facilities. Reference Materials should be stored on receipt as indicated on the label.

Real-Time Stability studies take place to ensure stated stability of this product

9. REFERENCES

10. ACKNOWLEDGEMENTS

11. FURTHER INFORMATION

Further information can be obtained as follows; This material: enquiries@nibsc.org WHO Biological Standards: http://www.who.int/biologicals/en/ JCTLM Higher order reference materials: http://www.bipm.org/en/committees/jc/jctlm/ Derivation of International Units: http://www.nibsc.org/standardisation/international_standards.aspx Ordering standards from NIBSC: http://www.nibsc.org/products/ordering.aspx NIBSC Terms & Conditions: http://www.nibsc.org/terms_and_conditions.aspx

12. CUSTOMER FEEDBACK

Customers are encouraged to provide feedback on the suitability or use of the material provided or other aspects of our service. Please send any comments to enquiries@nibsc.org

13. CITATION

In all publications, including data sheets, in which this material is referenced, it is important that the preparation's title, its status, the NIBSC code number, and the name and address of NIBSC are cited and cited correctly.

14. MATERIAL SAFETY SHEET

Classification in accordance with Directive 2000/54/EC, Regulation (EC) No 1272/2008: Not applicable or not classified

National Institute for Biological Standards and Control, Potters Bar, Hertfordshire, EN6 3QG. T +44 (0)1707 641000, nibsc.org WHO International Laboratory for Biological Standards, UK Official Medicines Control Laboratory

Medicines & Healthcare products Regulatory Agency



| Physical and Chemical properties | | | |
|---|--|----------------------------------|--------------------|
| Physical appearance: | | Corrosive: | No |
| Liquid | | | |
| Stable: Yes | | Oxidising: | No |
| Hygroscopi No | | Irritant: | No |
| с: | | | |
| Flammable: No | | Handling: See caution, Section 2 | |
| Other | | | |
| (specify): | | | |
| Toxicological properties | | | |
| roxideregidal propertied | | | |
| Effects of inhalation: Not | | established, a | void inhalation |
| Effects of ingestion: Not | | established, avoid ingestion | |
| Effects of skir | Not | established, | avoid contact with |
| absorption: | skin | | |
| Suggested First Aid | | | |
| ouggested first Ald | | | |
| Inhalation: See | Seek medical advice | | |
| Ingestion: See | Seek medical advice | | |
| Contact with Was | Wash with copious amounts of water. Seek | | |
| eyes: med | medical advice | | |
| Contact with Was | h thoro | oughly with wa | ter. |
| skin: | | | |
| Action on Spillage and Method of Disposal | | | |

Spillage of ampoule contents should be taken up with absorbent material wetted with an appropriate disinfectant. Rinse area with an appropriate disinfectant followed by water. Absorbent materials used to treat spillage should be treated as biological waste.

15. LIABILITY AND LOSS

In the event that this document is translated into another language, the English language version shall prevail in the event of any inconsistencies between the documents.

Unless expressly stated otherwise by NIBSC, NIBSC's Standard Terms and Conditions for the Supply of Materials (available at http://www.nibsc.org/About_Us/Terms_and_Conditions.aspx or upon request by the Recipient) ("Conditions") apply to the exclusion of all other terms and are hereby incorporated into this document by reference. The Recipient's attention is drawn in particular to the provisions of clause 11 of the Conditions.

16. INFORMATION FOR CUSTOMS USE ONLY

| Country of origin for customs purposes*: United Kingdom | | |
|--|--|--|
| * Defined as the country where the goods have been | | |
| produced and/or sufficiently processed to be classed as | | |
| originating from the country of supply, for example a change | | |
| of state such as freeze-drying. | | |
| Net weight: 4.5g | | |
| Toxicity Statement: Toxicity not assessed | | |
| Veterinary certificate or other statement if applicable. | | |
| Attached: No | | |