**CE Marked Material**  
**British Working Standard for Anti HIV-1**  
**NIBSC code: 99/750-xxx**  
**Instructions for use**  
**(Version 13.0, Dated 01/06/2020)**

**This material is an 'Annex II List A' IVD and complies with the requirements of the "EU in vitro diagnostic medical device directive 98/79/EC".**

**1. INTENDED USE**

**This product is CE marked for use as an IVD within the UK, EU member states and EEA countries. In all other territories this product can be used for research purposes only.**

The British Working Standard for Anti-HIV-1 (BWS) manufactured by NIBSC has been on the market since 1999. It serves as the UK Working Standard cited in the Guidelines for the Blood Transfusion Services in the United Kingdom, 7<sup>th</sup> Edition 2005.

The British Working Standard for Anti-HIV-1 (BWS) is intended for use in the field of *in vitro* diagnostics, in conjunction with diagnostic immunoassay test kits / systems for the detection of antibodies to HIV-1 to monitor the performance of these systems. It can be used to monitor the consistency of test performance using statistical process control on a daily basis and, over a period of time, as a retrospective monitor of batch performance. However it is for the user to establish suitability of purpose. It is expected that the BWS for Anti-HIV-1 will be detected by most 3<sup>rd</sup> generation assays. The British Working Standard for Anti-HIV-1 is for use in immunoassay test kits or systems that are capable of detecting antibodies to HIV-1

**2. CAUTION**

**This preparation is not for administration to humans or animals in the human food chain.**

This preparation contains material of human origin. It contains a dilution of 8 plasma samples known to be positive for antibodies to HIV-1. This plasma has been tested and found negative for HBsAg, and for antibodies to HCV, HIV-2 and HTLV-I / II, except for one sample that gave a very weak anti-HCV reaction when undiluted, although it was negative for HCV RNA by PCR. The final dilution of this sample in the BWS is 1 in 40,000.

As with all materials of biological origin, this preparation should be regarded as potentially hazardous to health. It should be used and discarded according to your own laboratory's safety procedures. Such safety procedures should include the wearing of protective gloves and avoiding the generation of aerosols. Care should be exercised in opening ampoules or vials, to avoid cuts.

**3. UNITAGE**

No unitage is assigned to this material.

**4. CONTENTS**

Country of origin of biological material: United Kingdom.  
Each vial of the British Working Standard for Anti-HIV-1 contains approximately 3ml of liquid and consists of:-

- Pooled and diluted anti-HIV-I-positive human plasma.
- Phosphate buffered saline (PBS-A).
- 5% Bovine serum albumin.
- 0.05% Bronidox as preservative.

**5. STORAGE**

Upon receipt, vials should be stored at +2-8°C. Reference materials are held at NIBSC within assured, temperature controlled storage facilities.

**6. DIRECTIONS FOR OPENING**

Vials have a screw cap; an internal stopper may also be present. The cap should be removed by turning anti-clockwise. Care should be taken to prevent loss of the contents. Please note: If a stopper is present on removal of the cap, the stopper should remain in the vial or be removed with the cap.

**7. USE OF MATERIAL**

**No attempt should be made to weigh out any portion of the material**

The user should follow the assay kit manufacturers' instructions. This preparation should be used in addition to the positive and negative controls supplied with each kit. Only the controls supplied with each kit by the manufacturer are used to determine the validity of assays and to calculate the cut-off value that forms the basis for donor screening or diagnosis.

**THIS PREPARATION IS SUPPLIED FOR USE IN ITS FINAL FORM AND MUST NOT BE FURTHER DILUTED OTHER THAN AS REQUIRED IN INDIVIDUAL TEST PROCEDURES.**

The Result Reporting System (RRS) has been developed by the National Institute for Biological Standards and Control (NIBSC) for the data monitoring of its serology and NAT quality control (QC) reagents. These include the Quality Control Reagent Unit (QCRU) and Clinical Virology Network (CVN) reagents. The system has been successfully running for serology assays for several years, collecting thousands of data points a year; and more recently for data derived from reagents used in Nucleic Acid-based Technologies (NAT) assays.  
<https://www.nibsc.org/products/rrs.aspx>

Users should be aware that different batches of the same assay kit and / or different assay conditions may give different values with this preparation.

**8. STABILITY**

The expiry date is stated on the vial label. Once opened, this preparation may be used for up to 3 months (as long as the expiry date is not exceeded), and should be stored at +2-8°C between each use.

The stability of this preparation is monitored by NIBSC. The stability of previous batches of the BWS have not lost potency when stored at +2-8°C for periods of up to 36 months. The shelf life assigned to this preparation is 36 months when stored at +2-8°C.

Users who have data supporting any deterioration in the characteristics of any Working Standard / Monitor Sample are encouraged to contact NIBSC.

**9. REFERENCES**

<http://www.transfusionguidelines.org.uk/index.asp?Publication=RB&Section=25&pageid=633>.

**10. ACKNOWLEDGEMENTS**

EC REP: Advena Ltd. Tower Business Centre, 2nd Flr., Tower Street, Swatar, BKR 4013 Malta

**11. FURTHER INFORMATION**

Further information can be obtained as follows;

This material: [enquiries@nibsc.org](mailto:enquiries@nibsc.org)

WHO Biological Standards:

<http://www.who.int/biologicals/en/>

JCTLM Higher order reference materials:

<http://www.bipm.org/en/committees/jc/jctlm/>

Derivation of International Units:

[http://www.nibsc.org/standardisation/international\\_standards.aspx](http://www.nibsc.org/standardisation/international_standards.aspx)

Ordering standards from NIBSC:

<http://www.nibsc.org/products/ordering.aspx>

NIBSC Terms & Conditions:

[http://www.nibsc.org/terms\\_and\\_conditions.aspx](http://www.nibsc.org/terms_and_conditions.aspx)

### 12. CUSTOMER FEEDBACK

Customers are encouraged to provide feedback on the suitability or use of the material provided or other aspects of our service. Please send any comments to [enquiries@nibsc.org](mailto:enquiries@nibsc.org)

### 13. CITATION

In all publications, including data sheets, in which this material is referenced, it is important that the preparation's title, its status, the NIBSC code number, and the name and address of NIBSC are cited and cited correctly.

### 14. MATERIAL SAFETY SHEET

Classification in accordance with Directive 2000/54/EC, Regulation (EC) No 1272/2008: Not applicable or not classified

Physical and Chemical properties	
Physical appearance: Liquid	Corrosive: No
Stable: Yes	Oxidising: No
Hygroscopic: No	Irritant: Yes (contains Bronidox)
Flammable: No	Handling: See caution, Section 2
Other (specify): Contains material of human origin.	
Toxicological properties	
Effects of inhalation: Not established. Avoid inhalation. May contain infectious HIV.	
Effects of ingestion: Not established. Avoid ingestion. May contain infectious HIV.	
Effects of skin absorption: Not established. Avoid contact with skin. May contain infectious HIV.	
Suggested First Aid	
Inhalation: Seek medical advice.	
Ingestion: Seek medical advice.	
Contact with eyes: Wash with copious amounts of water. Seek medical advice.	
Contact with skin: Wash thoroughly with water.	
Action on Spillage and Method of Disposal	
Spillage of vial contents should be taken up with absorbent material wetted with an appropriate disinfectant. Rinse area with an appropriate disinfectant followed by water. Absorbent materials used to treat spillage should be treated as biological waste.	

### 15. LIABILITY AND LOSS

In the event that this document is translated into another language, the English language version shall prevail in the event of any inconsistencies between the documents.

Unless expressly stated otherwise by NIBSC, NIBSC's Standard Terms and Conditions for the Supply of Materials (available at [www.nibsc.org/About\\_Us/Terms\\_and\\_Conditions.aspx](http://www.nibsc.org/About_Us/Terms_and_Conditions.aspx) or upon request by the Recipient) ("Conditions") apply to the exclusion of all other terms and are hereby incorporated into this document by reference. The Recipient's attention is drawn in particular to the provisions of clause

11 of the Conditions.

### 16. INFORMATION FOR CUSTOMS USE ONLY

**Country of origin for customs purposes\*:** United Kingdom

\* Defined as the country where the goods have been produced and/or sufficiently processed to be classed as originating from the country of supply, for example a change of state such as freeze-drying.

**Net weight:** 3g per vial

**Toxicity Statement:** Non-toxic.

**Veterinary certificate or other statement** if applicable.

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EC REP: Advena Ltd. Tower Business Centre, 2nd Flr., Tower Street, Swatara, BKR 4013 Malta

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[http://www.who.int/biologicals/reference\\_preparations/en/](http://www.who.int/biologicals/reference_preparations/en/)

Encargar patrones de NIBSC:

[http://www.nibsc.ac.uk/products/ordering\\_information/frequently\\_asked\\_questions.aspx](http://www.nibsc.ac.uk/products/ordering_information/frequently_asked_questions.aspx)

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[http://www.nibsc.ac.uk/terms\\_and\\_conditions.aspx](http://www.nibsc.ac.uk/terms_and_conditions.aspx)

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