



WHO International Standard
6th WHO International Standard for hepatitis C virus RNA for
nucleic acid amplification techniques
NIBSC code: 18/184
Instructions for use
(Version 3.0, Dated 27/11/2024)

5

1. INTENDED USE

The 6th WHO International Standard for hepatitis C virus (HCV) RNA for nucleic acid amplification techniques (NAT), NIBSC code 18/184, is intended to be used for the calibration of secondary reference reagents comprising HCV plasma used in HCV NAT [1]. The standard comprises lyophilized human plasma and HCV. The HCV was sourced from a high-titre HCV RNA genotype 1a-positive window-period donation. The standard has been lyophilized in 1.1 mL aliquots and stored at -20 °C. The material has been calibrated in International Units (IU), in parallel with the 5th WHO International Standard for HCV RNA for NAT, in a collaborative study involving 19 laboratories worldwide [2].

2. CAUTION

This preparation is not for administration to humans or animals in the human food chain.

None of the above

The preparation contains material of human origin. The HCV was sourced from a high-titre HCV RNA genotype 1a-positive window period donation. It was tested at the source laboratory and found to be negative for HCV antibody, HBsAg, HBV DNA, HIV-1/2 antibodies, and HIV RNA. It was also found to be negative for B19V DNA and HAV RNA in a pool of 256 donations. The pooled human plasma diluent was sourced from UK blood donations and had been tested and found negative for HIV antibody, HCV antibody, HBsAg and syphilis. It was also tested at NIBSC and found to be negative for B19V DNA and HAV RNA by NAT. As with all materials of biological origin, this preparation should be regarded as potentially hazardous to health. It should be used and discarded according to your own laboratory's safety procedures. Such safety procedures should include the wearing of protective gloves and avoiding the generation of aerosols. Care should be exercised in opening ampoules or vials, to avoid cuts.

3. UNITAGE

This material has been assigned a unitage of 257,000 IU/vial (5.41 log₁₀ IU/vial). Each vial must be reconstituted in 1.1 mL of nuclease-free water. Uncertainty: the assigned unitage does not carry an uncertainty associated with its calibration. The uncertainty may therefore be considered to be the variance of the vial weight after filling and was determined to be +/- 0.19%.

4. CONTENTS

Country of origin of biological material: United Kingdom.
Each vial contains 1.1 mL of lyophilized plasma containing infectious HCV.

5. STORAGE

Upon receipt, vials of lyophilized standard should be stored at -20 °C. This material is shipped on dry ice because stability assessment showed some instability upon storage at temperatures above +20 °C for 2 weeks. Please note because of the inherent stability of lyophilized material, NIBSC may ship these materials at ambient temperature.

6. DIRECTIONS FOR OPENING

Screw Cap Vials

7. USE OF MATERIAL

No attempt should be made to weigh out any portion of the freeze-dried material prior to reconstitution

The material should be reconstituted with 1.1 mL of deionized, nuclease-free molecular-grade water and left for a minimum of 20 minutes with occasional agitation before use. The reconstituted material has a final concentration of 233,636 IU/mL (~5.37 log₁₀ IU/mL). Once reconstituted, the International Standard should be diluted in the matrix appropriate to the assay e.g. in HCV RNA-negative plasma, and should be extracted prior to HCV RNA measurement. The International Standard should be used to calibrate secondary reference materials, for example, by determining the equivalent concentration of secondary reference reagent being calibrated, against the International Standard, in parallel, as described elsewhere [1]. The secondary reference reagent can then be assigned a concentration in IU.

8. STABILITY (Add or amend as necessary)

Reference materials are held at NIBSC within assured, temperature-controlled storage facilities. Reference Materials should be stored on receipt as indicated on the label.

It is the policy of WHO not to assign an expiry date to their international reference materials. Accelerated degradation studies have indicated that this material is suitably stable, when stored at -20 °C, for the assigned values to remain valid until the material is withdrawn or replaced. Users who have data supporting any deterioration in the characteristics of any reference preparation are encouraged to contact NIBSC.

NIBSC follows the policy of WHO with respect to its reference materials.

9. REFERENCES

1. WHO manual for the preparation of secondary reference materials for in vitro diagnostic assays designed for infectious disease nucleic acid or antigen detection: calibration to WHO International Standards. WHO Technical Report Series 2017. Geneva, Switzerland:WHO 2017; 1004,389-455.
2. Fryer JF, Rigsby P, Hockley JG, Morris CL, and the Collaborative Study Group. A Collaborative Study to Evaluate the Proposed 6th WHO International Standard for Hepatitis C Virus (HCV) RNA for Nucleic Acid Amplification Techniques (NAT). WHO ECBS Report 2019; WHO/BS/2019.2358.

10. ACKNOWLEDGEMENTS

We gratefully acknowledge the important contributions of the collaborative study participants.

11. FURTHER INFORMATION

Further information can be obtained as follows;
This material: enquiries@nibsc.org
WHO Biological Standards:

<http://www.who.int/biologicals/en/>
 JCTLM Higher order reference materials:
<http://www.bipm.org/en/committees/jc/jctlm/>
 Derivation of International Units:
http://www.nibsc.org/standardisation/international_standards.aspx
 Ordering standards from NIBSC:
<http://www.nibsc.org/products/ordering.aspx>
 NIBSC Terms & Conditions:
http://www.nibsc.org/terms_and_conditions.aspx

12. CUSTOMER FEEDBACK

Customers are encouraged to provide feedback on the suitability or use of the material provided or other aspects of our service. Please send any comments to enquiries@nibsc.org

13. CITATION

In all publications, including data sheets, in which this material is referenced, it is important that the preparation's title, its status, the NIBSC code number, and the name and address of NIBSC are cited and cited correctly.

14. MATERIAL SAFETY SHEET (Add or amend as necessary)

Classification in accordance with Directive 2000/54/EC, Regulation (EC) No 1272/2008: Not applicable or not classified

| Physical and Chemical properties | |
|---|---|
| Physical appearance: Lyophilized powder | Corrosive: No |
| Stable: Yes | Oxidising: No |
| Hygroscopic: No | Irritant: No |
| Flammable: No | Handling: See caution, Section 2 |
| Other (specify): | Contains infectious hepatitis C virus and human plasma |
| Toxicological properties | |
| Effects of inhalation: | Avoid. Contains infectious hepatitis C virus and human plasma |
| Effects of ingestion: | Avoid. Contains infectious hepatitis C virus and human plasma |
| Effects of skin absorption: | Avoid. Contains infectious hepatitis C virus and human plasma |
| Suggested First Aid | |
| Inhalation: | Seek medical advice |
| Ingestion: | Seek medical advice |
| Contact with eyes: | Wash with copious amounts of water. Seek medical advice |
| Contact with skin: | Wash thoroughly with water. |
| Action on Spillage and Method of Disposal | |
| Spillage of ampoule contents should be taken up with absorbent material wetted with an appropriate disinfectant. Rinse area with an appropriate disinfectant followed by water. Absorbent materials used to treat spillage should be treated as biological waste. | |

15. LIABILITY AND LOSS

In the event that this document is translated into another language, the English language version shall prevail in the event of any inconsistencies between the documents.

Unless expressly stated otherwise by NIBSC, NIBSC's Standard Terms and Conditions for the Supply of Materials (available at http://www.nibsc.org/About_Us/Terms_and_Conditions.aspx or upon request by the Recipient) ("Conditions") apply to the exclusion of all other terms and are hereby incorporated into this document by reference. The Recipient's attention is drawn in particular to the provisions of clause 11 of the Conditions.

16. INFORMATION FOR CUSTOMS USE ONLY

| |
|---|
| Country of origin for customs purposes*: United Kingdom * Defined as the country where the goods have been produced and/or sufficiently processed to be classed as originating from the country of supply, for example a change of state such as freeze-drying. |
| Net weight: 1.1 g |
| Toxicity Statement: Non-toxic |
| Veterinary certificate or other statement if applicable. |
| Attached: No |

17. CERTIFICATE OF ANALYSIS

NIBSC does not provide a Certificate of Analysis for WHO Biological Reference Materials because they are internationally recognised primary reference materials fully described in the instructions for use. The reference materials are established according to the WHO Recommendations for the preparation, characterization and establishment of international and other biological reference standards [https://www.who.int/publications/m/item/annex2-trs932\(revised 2004\)](https://www.who.int/publications/m/item/annex2-trs932(revised%202004)). They are officially endorsed by the WHO Expert Committee on Biological Standardization (ECBS) based on the report of the international collaborative study which established their suitability for the intended use.