Medicines & Healthcare products Regulatory Agency



#### Working Standard Anti-D Immunoglobulin Control Preparation NIBSC code: 16/278 Instructions for use (Version 1.0, Dated 26/03/2019)

#### This material is not for in vitro diagnostic use.

#### 1. INTENDED USE

Preparation 16/278 is intended for use as a control in anti-D potency assays of immunoglobulin products to monitor assay consistency. It is not intended as a standard or reference preparation. For this purpose, the International Standard for anti-D immunoglobulin, 16/332 should be used.

#### 2. CAUTION

# This preparation is not for administration to humans or animals in the human food chain.

The preparation contains material of human origin, and either the final product or the source materials, from which it is derived, have been tested and found negative for HBsAg, anti-HIV and HCV RNA. As with all materials of biological origin, this preparation should be regarded as potentially hazardous to health. It should be used and discarded according to your own laboratory's safety procedures. Such safety procedures should include the wearing of protective gloves and avoiding the generation of aerosols. Care should be exercised in opening ampoules or vials, to avoid cuts.

# 3. UNITAGE

N/A

#### 4. CONTENTS

Country of origin of biological material: USA. Clinical grade anti-D immunoglobulin was distributed into ampoules and lyophilized.

#### 5. STORAGE

Store unopened ampoules at -20°C or below. Please note: because of the inherent stability of lyophilized material, NIBSC may ship these materials at ambient temperature.

## 6. DIRECTIONS FOR OPENING

DIN ampoules have an 'easy-open' coloured stress point, where the narrow ampoule stem joins the wider ampoule body. Various types of ampoule breaker are available commercially. To open the ampoule, tap the ampoule gently to collect material at the bottom (labelled) end and follow manufactures instructions provided with the ampoule breaker.

#### 7. USE OF MATERIAL

No attempt should be made to weigh out any portion of the freeze-dried material prior to reconstitution

Reconstitute the ampoule contents with 1.0 ml of distilled water containing 0.02% (w/v) sodium azide.

ALLOW AT LEAST 10 MIN FOR RECONSTITUTION.

The total protein concentration of the reconstituted contents is approximately 50mg/ml. The intended use is as an assay run control using the same dilution series as the 3<sup>rd</sup> IS for anti-D immunoglobulin (16/332). End users will need to assign a mean value and acceptance limits for their laboratory method.

#### 8. STABILITY

NIBSC follows the policy of WHO with respect to its reference materials. It is the policy of WHO not to assign an expiry date to its reference materials. They remain valid with the assigned potency and status until withdrawn or amended.

Reference materials are held at NIBSC within assured, temperaturecontrolled storage facilities. Reference Materials should be stored on receipt as indicated on the label.

Once reconstituted, diluted or aliquoted, users should determine the stability of the material according to their own method of preparation, storage and use.

Users who have data supporting any deterioration in the characteristics of any reference preparation are encouraged to contact NIBSC. For information specific to a particular biological standard, contact NIBSC.

#### 9. REFERENCES

#### **10. ACKNOWLEDGEMENTS**

#### 11. FURTHER INFORMATION

Further information can be obtained as follows; This material: enquiries@nibsc.org WHO Biological Standards: http://www.who.int/biologicals/en/ JCTLM Higher order reference materials: http://www.bipm.org/en/committees/jc/jctlm/ Derivation of International Units: http://www.nibsc.org/standardisation/international\_standards.aspx Ordering standards from NIBSC: http://www.nibsc.org/products/ordering.aspx NIBSC Terms & Conditions: http://www.nibsc.org/terms\_and\_conditions.aspx

## 12. CUSTOMER FEEDBACK

Customers are encouraged to provide feedback on the suitability or use of the material provided or other aspects of our service. Please send any comments to enquiries@nibsc.org

#### 13. CITATION

In all publications, including data sheets, in which this material is referenced, it is important that the preparation's title, its status, the NIBSC code number, and the name and address of NIBSC are cited and cited correctly.

#### 14. MATERIAL SAFETY SHEET

Classification in accordance with Directive 2000/54/EC, Regulation (EC) No 1272/2008: Not applicable or not classified

Physical and Chemical properties				
Physical appearance:			Corrosive:	No
Lyophilisate				
Stable:	Yes		Oxidising:	No
Hygroscopic:	No		Irritant:	No
Flammable:	No		Handling: Se	e caution, Section 2
Other (specify): See caution, section 2				
Toxicological properties				
Effects of inhalation:		Not established, avoid inhalation		
Effects of ingestion:		Not established, avoid ingestion		
Effects of skin absorption:		Not established, avoid contact with skin		
Suggested First Aid				
Inhalation: Seek medical advice				
Ingestion: Seek medical advice				

National Institute for Biological Standards and Control,

Potters Bar, Hertfordshire, EN6 3QG. T +44 (0)1707 641000, nibsc.org WHO International Laboratory for Biological Standards,

UK Official Medicines Control Laboratory

# Medicines & Healthcare products Regulatory Agency



Contact with skin:	Mash thoroughly with water.
Contact with eyes:	Wash with copious amounts of water. Seek medical advice

#### Action on Spillage and Method of Disposal

Spillage of contents should be taken up with absorbent material wetted with an appropriate disinfectant. Rinse area with an appropriate disinfectant followed by water.

Absorbent materials used to treat spillage should be treated as biological waste.

# 15. LIABILITY AND LOSS

In the event that this document is translated into another language, the English language version shall prevail in the event of any inconsistencies between the documents.

Unless expressly stated otherwise by NIBSC, NIBSC's Standard Terms and Conditions for the Supply of Materials (available at http://www.nibsc.org/About\_Us/Terms\_and\_Conditions.aspx or upon request by the Recipient) ("Conditions") apply to the exclusion of all other terms and are hereby incorporated into this document by reference. The Recipient's attention is drawn in particular to the provisions of clause 11 of the Conditions.

# 16. INFORMATION FOR CUSTOMS USE ONLY

Country of origin for customs purposes\*: United Kingdom \* Defined as the country where the goods have been produced and/or sufficiently processed to be classed as originating from the country of supply, for example a change of state such as freeze-drying. Net weight: 0.1g Toxicity Statement: Toxicity not assessed

Veterinary certificate or other statement if applicable. Attached: No

National Institute for Biological Standards and Control, Potters Bar, Hertfordshire, EN6 3QG. T +44 (0)1707 641000, nibsc.org WHO International Laboratory for Biological Standards, UK Official Medicines Control Laboratory