Influenza Reagent
Influenza Virus Infectious
A/Indiana/02/2020 (H1N1)
NIBSC code: 22/138
Instructions for use
(Version 1.0, Dated 18/08/2022)

§

1. INTENDED USE

Reagent 22/138 is prepared from A/Indiana/02/2020 (H1N1), which was processed in 250µl volumes as liquid stock. The known passage history of A/Indiana/02/2020 is attached.

2. CAUTION

The material is not of human or bovine origin. This preparation is not for administration to humans or animals

As with all materials of biological origin, this preparation should be regarded as potentially hazardous to health. It should be used and discarded according to your own laboratory's safety procedures. Such safety procedures should include the wearing of protective gloves and avoiding the generation of aerosols. Care should be exercised in opening ampoules or vials, to avoid cuts.

3. UNITAGE

No unitage is assigned to this material.

4. CONTENTS

Country of origin of biological material: United Kingdom. Each vial contains 250µl (nominal) of infectious influenza virus as allantoic fluid from SPF embryonated hen's eggs.

5. STORAGE

Store in the dark at -70°C or below.

Material type: Liquid – will be shipped according to the storage and shipping conditions of the product

6. DIRECTIONS FOR OPENING

Vials have a screw cap; an internal stopper may also be present. The cap should be removed by turning anti-clockwise. Care should be taken to prevent loss of the contents. Please note: If a stopper is present on removal of the cap, the stopper should remain in the vial or be removed with the cap.

7. USE OF MATERIAL

Ready to use.

8. STABILITY

Reference materials are held at NIBSC within assured, temperature-controlled storage facilities. Reference Materials should be stored on receipt as indicated on the label.

NIBSC follows the policy of WHO with respect to its reference materials.

9. REFERENCES

NA

10. ACKNOWLEDGEMENTS

NA



11. FURTHER INFORMATION

Further information can be obtained as follows;

This material: enquiries@nibsc.org WHO Biological Standards:

http://www.who.int/biologicals/en/

JCTLM Higher order reference materials:

http://www.bipm.org/en/committees/jc/jctlm/

Derivation of International Units:

 $http://www.nibsc.org/standardisation/international_standards.aspx$

Ordering standards from NIBSC:

http://www.nibsc.org/products/ordering.aspx

NIBSC Terms & Conditions:

http://www.nibsc.org/terms_and_conditions.aspx

12. CUSTOMER FEEDBACK

Customers are encouraged to provide feedback on the suitability or use of the material provided or other aspects of our service. Please send any comments to enquiries@nibsc.org

13. CITATION

In all publications, including data sheets, in which this material is referenced, it is important that the preparation's title, its status, the NIBSC code number, and the name and address of NIBSC are cited and cited correctly.

14. MATERIAL SAFETY SHEET

Classification in accordance with Directive 2000/54/EC, Regulation (EC) No 1272/2008: Not applicable or not classified

(EC) No 1272/2008: Not	t applicable or not classified				
Physical and Chemical properties					
Physical appearance: Clear liquid	Corrosive: No				
Stable: Yes	Oxidising: No				
Hygroscopi No c:	Irritant: No				
Flammable: No	Handling: See caution, Section 2				
Other Live influenza virus (specify):					
Toxicological properties					
Effects of inhalation:	Likelihood of influenza virus infection.				
Effects of ingestion:	Not established, avoid ingestion				
Effects of skin	Not established, avoid contact with				
absorption:	skin				
Suggested First Aid					
Inhalation: Seek medical advice					
Ingestion: Seek medical advice					
Contact with Wasl	with copious amounts of water. Seek				
eyes: med	al advice				
Contact with Wash thoroughly with water. kin:					
Action on Spillage and Method of Disposal					

Spillage of ampoule contents should be taken up with absorbent material wetted with an appropriate disinfectant.

Rinse area with an appropriate disinfectant followed by water. Absorbent materials used to treat spillage should be treated as biological waste.

National Institute for Biological Standards and Control, Potters Bar, Hertfordshire, EN6 3QG. T +44 (0)1707 641000, nibsc.org WHO International Laboratory for Biological Standards, UK Official Medicines Control Laboratory



15. LIABILITY AND LOSS

In the event that this document is translated into another language, the English language version shall prevail in the event of any inconsistencies between the documents.

Unless expressly stated otherwise by NIBSC, NIBSC's Standard Terms and Conditions for the Supply of Materials (available at http://www.nibsc.org/About_Us/Terms_and_Conditions.aspx or upon request by the Recipient) ("Conditions") apply to the exclusion of all other terms and are hereby incorporated into this document by reference. The Recipient's attention is drawn in particular to the provisions of clause 11 of the Conditions.

16. INFORMATION FOR CUSTOMS USE ONLY

Country of origin for customs purposes*: United Kingdom
* Defined as the country where the goods have been
produced and/or sufficiently processed to be classed as
originating from the country of supply, for example a change
of state such as freeze-drying.

Net weight: 0.25g per vial.

Toxicity Statement: Non-toxic

Veterinary certificate or other statement if applicable.

Attached: No

Passage history of A/Indiana/02/2020 (H1N1)

Cumulative number of passages	Passage numbers at each stage	Lot	Laboratory
E3	E3	3026055186	CDC, USA
E5	E3/E2	46950	NIBSC, UK

Sterility: No visible contamination was detected in a variety of media (tryptose soya broth, thioglycolate broth, Sabouraud's broth and blood agar plates) after 14 days incubation.

The HA and NA sequence of this virus is available at GISAID with the accession number EPI_ISL_13453558.