



**Non WHO Reference Material
Rec. Staphylokinase (STAR) NIBSC Reagent
NIBSC code: 94/718
Instructions for use
(Version 6.0, Dated 04/01/2011)**

This material is not for in vitro diagnostic use.

1. INTENDED USE

This preparation is intended for use as a standard in plasminogen activation assays designed to measure the potency of staphylokinase solutions

2. CAUTION

This preparation is not for administration to humans or animals in the human food chain.

The preparation contains material of human origin, and either the final product or the source materials, from which it is derived, have been tested and found negative for HBsAg, anti-HIV and HCV RNA. As with all materials of biological origin, this preparation should be regarded as potentially hazardous to health. It should be used and discarded according to your own laboratory's safety procedures. Such safety procedures should include the wearing of protective gloves and avoiding the generation of aerosols. Care should be exercised in opening ampoules or vials, to avoid cuts.

3. UNITAGE

Each ampoule contains 1.0 unit of recombinant Staphylokinase (STAR). The assigned content was valid at the time of manufacture-no information available on long term stability.

4. CONTENTS

Country of origin of biological material: United Kingdom.
The source material was supplied as a frozen solution of 25ml of STAR concentrate (1.05mg/ml) in a solution of 15mM potassium phosphate, 150mM NaCl, pH 7.4. This frozen solution was thawed and diluted with a 0.01M phosphate buffer (pH 7.4) containing 5mgs of human albumin per ml to a final concentration of 15 micrograms of STAR per ml. This batch of ampoules, coded 94/718, contains the residue after freeze drying of an accurately measured 1ml of a Staphylokinase solution (15 micrograms STAR in 10mM phosphate buffer, pH7.4 containing 5mg human albumin).

5. STORAGE

Unopened ampoules should be stored in the dark at or below -20°C.
Please note: because of the inherent stability of lyophilized material, NIBSC may ship these materials at ambient temperature.

6. DIRECTIONS FOR OPENING

DIN ampoules have an 'easy-open' coloured stress point, where the narrow ampoule stem joins the wider ampoule body.
Tap the ampoule gently to collect the material at the bottom (labeled) end. Ensure that the disposable ampoule safety breaker provided is pushed down on the stem of the ampoule and against the shoulder of the ampoule body. Hold the body of the ampoule in one hand and the disposable ampoule breaker covering the ampoule stem between the thumb and first finger of the other hand. Apply a bending force to open the ampoule at the coloured stress point, primarily using the hand holding the plastic collar.
Care should be taken to avoid cuts and projectile glass fragments that might enter the eyes, for example, by the use of suitable gloves and an eye shield. Take care that no material is lost from the ampoule and no glass falls into the ampoule. Within the ampoule is dry nitrogen gas at slightly less than atmospheric pressure. A new disposable ampoule breaker is provided with each DIN ampoule.

Care should be taken to avoid cuts and projectile glass fragments that might enter the eyes, for example, by the use of suitable gloves and an eye shield. Take care that no material is lost from the ampoule and no glass falls into the ampoule. Within the ampoule is dry nitrogen gas at slightly less than atmospheric pressure. A new disposable ampoule breaker is provided with each DIN ampoule.

7. USE OF MATERIAL

No attempt should be made to weigh out any portion of the freeze-dried material prior to reconstitution

The total contents of the ampoule should be reconstituted with 1.0ml distilled water and dissolved gently by swirling to avoid froth. The reconstituted reagent should be used as soon as possible after reconstitution

8. STABILITY

Reference materials are held at NIBSC within assured, temperature-controlled storage facilities. Reference Materials should be stored on receipt as indicated on the label.

NIBSC follows the policy of WHO with respect to its reference materials.

Following ampouling by procedures outlined by Campbell (1974) the stability was established by assaying the contents of ampoules stored at elevated temperatures (Jerne & Perry, 1956). It was observed that the reagent was stable for 18 months at 37°C using a modified clot lysis assay (Heath & Gaffney, 1990).

9. REFERENCES

Campbell, P.J., (1974) J. Biol. Stand. 2: 259-267

Heath, A. Gaffney, P.J, (1990) Thromb. Haemost. 64, 267-269

Jerne, N.K, Perry, W.L.M. (1956) Bull. WHO 14: 167-182.

10. ACKNOWLEDGEMENTS

We are grateful to the Hans Knoll Institute and the good offices of Professor D. Collen, for supplying the staphylokinase source material.

11. FURTHER INFORMATION

Further information can be obtained as follows;
This material: enquiries@nibsc.org
WHO Biological Standards:
<http://www.who.int/biologicals/en/>
JCTLM Higher order reference materials:
<http://www.bipm.org/en/committees/jc/jctlm/>
Derivation of International Units:
http://www.nibsc.org/standardisation/international_standards.aspx
Ordering standards from NIBSC:
<http://www.nibsc.org/products/ordering.aspx>
NIBSC Terms & Conditions:
http://www.nibsc.org/terms_and_conditions.aspx

12. CUSTOMER FEEDBACK

Customers are encouraged to provide feedback on the suitability or use of the material provided or other aspects of our service. Please send any comments to enquiries@nibsc.org

13. CITATION

In all publications, including data sheets, in which this material is referenced, it is important that the preparation's title, its status, the NIBSC code number, and the name and address of NIBSC are cited and cited correctly.

14. MATERIAL SAFETY SHEET

Classification in accordance with Directive 2000/54/EC, Regulation (EC)
No 1272/2008: Not applicable or not classified

Physical and Chemical properties



Physical appearance: Freeze-dried powder	Corrosive: No
Stable: Yes	Oxidising: No
Hygroscopic: Yes	Irritant: No
Flammable: No	Handling: See caution, Section 2
Other (specify): Contains material of human origin	
Toxicological properties	
Effects of inhalation: Not established, avoid inhalation	
Effects of ingestion: Not established, avoid ingestion	
Effects of skin absorption: Not established, avoid contact with skin	
Suggested First Aid	
Inhalation: Seek medical advice	
Ingestion: Seek medical advice	
Contact with eyes: Wash with copious amounts of water. Seek medical advice	
Contact with skin: Wash thoroughly with water.	
Action on Spillage and Method of Disposal	
Spillage of ampoule contents should be taken up with absorbent material wetted with an appropriate disinfectant. Rinse area with an appropriate disinfectant followed by water. Absorbent materials used to treat spillage should be treated as biological waste.	

15. LIABILITY AND LOSS

In the event that this document is translated into another language, the English language version shall prevail in the event of any inconsistencies between the documents.

Unless expressly stated otherwise by NIBSC, NIBSC's Standard Terms and Conditions for the Supply of Materials (available at http://www.nibsc.org/About_Us/Terms_and_Conditions.aspx or upon request by the Recipient) ("Conditions") apply to the exclusion of all other terms and are hereby incorporated into this document by reference. The Recipient's attention is drawn in particular to the provisions of clause 11 of the Conditions.

16. INFORMATION FOR CUSTOMS USE ONLY

Country of origin for customs purposes*: United Kingdom * Defined as the country where the goods have been produced and/or sufficiently processed to be classed as originating from the country of supply, for example a change of state such as freeze-drying.
Net weight: 10mg
Toxicity Statement: Toxicity not assessed
Veterinary certificate or other statement if applicable. Attached: No