

## **Centre for AIDS Reagents.**



## **Data Sheet**

NAME:	U87.CD4 cell lines (Chemokine receptor transfectants).
MANIE.	667.664 cen mes (enemokine receptor transfectants)

**REPOSITORY REFERENCE:** ARP069-073

**SPECIES/TYPE:** Human glioma cells expressing human chemokine receptors..

**DESCRIPTION:** Adherent cells expressing chemokine co-receptors on cell surface.

ARP:	Chemokine Co-receptor
069 Control	U87.CD4
069	CCR1
070	CCR2b
071	CCR3
072	CCR5
073	CXCR4

CD4 was introduced with pMV7neo retroviral vector conferring G418 resistance. All co-receptor genes were introduced with retroviral vector pBABEpuro.U87-CXCR4 has no detectable cell surface CD4 expression. All cells except 069 Control are positive for Amphotrophic Murine Leukaemia Virus. All cells are negative for Mycoplasma.

**SPECIAL CHARACTERISTICS:** These cells can be used for HIV/SIV infections for virus time

course/production assays. All U87 cells endogenously express the GPR1 orphan receptor and ARP071-073 cells also express the Bonzo/STRL33 receptor. These receptors can be used,

primarily by certain SIV isolates for entry.

**GROWTH CHARACTERISTICS:** Cells divide slowly; however, they should not be allowed to

become over-confluent. In general, pass cells at 80% confluence and split no more than 1:10. All cells are sensitive

to acidic medium.

Version 1 Page 1 of 2

**CULTURE MEDIUM: DMEM**, 85% Heat inactivated Foetal calf serum, 15% G418, 300µg/ml Puromycin 1µg/ml. (Do not add Puromycin to 069c). Supplement with penicillin and streptomycin. STORAGE: Liquid nitrogen **SOURCE:** Dr Littman and Dr HongKui Deng REFERENCES: Bjorndal A et al, 1997, J Virol 71: 7478-7487 **ACKNOWLEDGEMENTS:** Publications should acknowledge the donor of the reagent and the Programme EVA Centre for AIDS Reagents. Suggested wording can be found on our website http://www.nibsc.ac.uk/spotlight/aidsreagent/index.html in the "Acknowledgements" section. Please also ensure that you send us a copy of any papers resulting from work using reagents acquired through CFAR (this can be electronically or as a paper copy)

NOTE:

Enquiries from commercial organisations should contact the New York University Industrial Liaison Office (Tel: +212 263

8178) to negotiate licensing agreements.

Version 1 Page 2 of 2