

**Centre for AIDS Reagents.** 

## Data Sheet.



NAME :	HIV-1 MC99IIIB∆Tat-Rev
<b>REPOSITORY REFERENCE :</b>	ARP170
ORIGINAL SOURCE :	HIV-1 was generated by transfection of the proviral plasmid IIIB, an HXB-3 derivative. Oligonucleotide-directed mutagenesis with the bacteriophage M13 system (Amersham) was used to introduce five point mutations into the coding regions of <i>tat</i> and <i>rev</i> .
SPECIAL CHARACTERISTICS :	This virus will only propagate in CEM-TART and IA2 cells. Mutant virus produced using these cells allows studies of the HIV-1 life cycle in a biologically contained system.
PROPAGATION :	CEM-TART and IA2 cells (ARP056)
STERILITY :	Previously tested positive for Mycoplasma. Treated with Ciprofloxacin, currently awaiting results to confirm negative status.
STORAGE :	Liquid Nitrogen.
SOURCE :	Drs Herbert Chen, Terence Boyle, Michael Malim, Brian Cullen and H Kim Lyerly. (courtesy of the NIH AIDS Research and Reference Reagent Program).

**REFERENCE :** 

Chen H et al (1992), Proc. Natl. Acad. Sci USA <u>89</u>:7678-7682.

## ACKNOWLEDGEMENTS :

Publications should acknowledge the donor of the reagent and the Programme EVA Centre for AIDS Reagents. Suggested wording can be found on our website at <u>http://www.nibsc.ac.uk/spotlight/aidsreagent/index.html</u> in the "Acknowledgements" section.

Please also ensure that you send us a copy of any papers resulting from work using reagents acquired through CFAR (this can be electronically or as a paper copy)