

Data Sheet

DESCRIPTION :	HIV-1 Neutralizing Serum #1
REPOSITORY REFERENCE :	ARP509
HOST :	Human
TITRE :	Strain dependent
CHARACTERISTICS :	<p>Sample is pooled serum of three separate bleeds from the same patient. An international collaborative study between six laboratories each using their own neutralisation assays demonstrated cross-reactivity with other HIV-1 strains and the mean neutralisation titres for some strains were:</p> <p>HIV IIB/LAV = 1:49; HIV RF = 1:63; HIV MN = 1:804; HIV SF2 = 1:329.</p>
PRESENTATION :	1 ampoule lyophilized protein. Reconstitute with 0.25 ml water.
SOURCE :	Dr L Vujcic
ACKNOWLEDGEMENTS :	<p>Publications should acknowledge the donor of the reagent and the Programme EVA Centre for AIDS Reagents. Suggested wording can be found on our website at http://www.nibsc.ac.uk/spotlight/aidsreagent/index.html in the "Acknowledgements" section.</p> <p>Please also ensure that you send us a copy of any papers resulting from work using reagents acquired through CFAR (this can be electronically or as a paper copy)</p>