



Centre for AIDS Reagents

Data Sheet

| NAME: | SIV gag p27 (aa233-247) |
|------------------------------|---|
| REPOSITORY REFERENCE: | ARP706 |
| SEQUENCE: | GSDIAGTTSSVDEQI |
| VIRUS STRAIN: | SIVmac251 |
| REGION OF HOMOLOGY: | Gag p27 (233-247 Los Alamos database) |
| SIZE: | 15mer |
| SPECIAL CHARACTERISTICS: | Possible T-cell epitope |
| SOLUBILITY: | Soluble in dilute base (i.e. 0.1% NH ₄ OH) or may require DMSO to dissolve, then dilute out with water |
| ACKNOWLEDGEMENTS: | Publications should acknowledge the donor of the reagent and the Programme EVA Centre for AIDS Reagents. Suggested wording can be found on our website in the "Acknowledgement" section at:- |
| | www.nibsc.ac.uk/spotlight/centre_for_aids_reagents.aspx |
| | Please also ensure that you send us a copy of any papers resulting from work using reagents acquired through CFAR, this can be by e-mail or printed copy. |

