

Centre for AIDS Reagents

Data Sheet

NAME: HIV-1 env gp120 V3 (GUN)

REPOSITORY REFERENCE: **ARP770-771**

VIRUS STRAIN: HIV-1 GUN

REGION OF HOMOLOGY: env gp120 V3 loop (305-318)

SIZE: 14mers

SPECIAL CHARACTERISTICS: Two V3 loop variants of Japanese HIV-1 GUN

SOLUBILITY: Readily soluble in water

ACKNOWLEDGEMENTS : Publications should acknowledge the donor of the reagent and the Programme EVA Centre for AIDS Reagents. Suggested wording can be found on our website at <http://www.nibsc.ac.uk/spotlight/aidsreagent/index.html> in the "Acknowledgements" section. Please also ensure that you send us a copy of any papers resulting from work using reagents acquired through CFAR (this can be electronically or as a paper copy)

SEQUENCE :

ARP770	KSITIGPGRAFHAI (Parent Strain)
ARP771	KSITIGSGRAFHAI (Variant Strain)