

Centre for AIDS Reagents

Data Sheet

NAME:	C-HIV-2 KARPAS
REPOSITORY REFERENCE:	ARP1167
PRESENTATION:	1 ml Virus infected supernatant
CELL LINE FOR PROPAGATION:	KARPAS cells
CULTURE MEDIUM:	RPMI 1640 + 10% FCS
STORAGE:	-80°C
SOURCE:	Dr. Abraham Karpas
REFERENCE:	Lancet. 1987 Jul 18;2(8551):132-3.
ACKNOWLEDGEMENTS:	<p>Publications should acknowledge the donor of the reagent and the Programme EVA Centre for AIDS Reagents. Suggested wording can be found on our website in the “Acknowledgements” section at:</p> <p>www.nibsc.ac.uk/spotlight/centre_for_aids_reagents.aspx</p> <p>Please also ensure that you send us a copy of any papers resulting from work using reagents acquired through CFAR (this can be electronically or as a paper copy)</p>