

Data Sheet

NAME:	pTRJO4551 clone 58 (SVPB17)
REPOSITORY REFERENCE	ARP2060
NOTE:	This clone is also available as a member of a panel set, see ARP2066
PROVIDED:	20 µg plasmid DNA/per vial (0.47 mg/ml)
HOST STRAIN:	MAX Efficiency STBL2™
CLONING SITE:	The env/rev cassette was TA cloned into pcDNA3.1/V5-His© TOPO®, and colonies were screened to obtain env/rev clones in the correct orientation with the CMV promoter. The size of the insert is 2982 bp.
CLONING VECTOR:	pcDNA3.1/V5-His© TOPO®. The size of the cloning vector including the insert is 8505 bp.
DESCRIPTION:	A PCR fragment containing full-length env and rev genes was derived from plasma virion-associated RNA from a subject acutely infected with a clade B virus by reverse transcription and nested PCR amplification procedures. The env/rev cassette was cloned into the pcDNA3.1/V5-His© TOPO® expression vector by the direct insertion of PCR products into a plasmid vector. A single transformed ampicillin-resistant <i>E. coli</i> colony was selected and the insert and correct orientation with CMV promoter was proved by the sequencing.
SPECIAL CHARACTERISTICS:	The clone represents env/rev cassette from a subject with acute subtype B infection (male to male transmission). The clone expresses a functional env/rev cassette and can be used to generate pseudotyped infectious virions. pTRJO4551.58 Env containing pseudovirions are included in a standard virus neutralization panel for subtype B strains (SVPB17).

PLASMID EXPANSION:

It is recommended that this plasmid be expanded using MAX Efficiency STBL2™ Competent Cells in LB medium at 30°C.

GENE BANK:

Accession number is AY835450

STORAGE:

-80°C

SOURCE:

Drs. B. H. Hahn, X. Wei and G. M. Shaw (Courtesy of NIH AIDS Research and Reference Reagent Programme.)

REFERENCE :

Li, M., Gao F., Mascola J.R., Stamatatos L., Polonis V.R., Koutsoukos M., Voss G., Goepfert P., Gilbert P., Greene K.M., Bilska M., Kothe D.L., Salazar- Gonzalez J.F., Wei X., Decker J.M., Hahn B.H., and Montefiori D.C. Human immunodeficiency virus type 1 env clones from acute and early subtype B infections for standardized assessments of vaccine-elicited neutralizing antibodies. *J. Virology* **79**(16): 10108-10125, 2005.

ACKNOWLEDGEMENTS:

Publications should acknowledge the donor of the reagent and the Programme EVA Centre for AIDS Reagents. Suggested wording can be found on our website at <http://www.nibsc.ac.uk/spotlight/aidsreagent/index.html> in the “Acknowledgements” section.

Please also ensure that you send us a copy of any papers resulting from work using reagents acquired through CFAR (this can be electronically or as a paper copy)

NOTE:

Scientists at for-profit institutions or who intend commercial use of this reagent must contact William S. White, UAB Research Foundation, The Office of Intellectual Property Management, AB 1120G, 1530 3rd Ave. S, Birmingham AL 35294-0111, Tel: 205-996-2550 Fax: 205-934-5427, email: wswhite@uab.edu, before the reagent can be released.